Cover page

Arabic..... 2

Chinese 3

Czech 4

Dutch 5

English 6

Finnish 7

French..... 8

German..... 9

Greek 10

Italian..... 11

Polish 12

Romanian 13

Russian..... 14

Slovenian 15

Spanish 16

Turkish 17

Ukrainian 18

Portuguese 19

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Buig de sensor niet.

English

3-Lead ECG Monitoring Electrode Array

The A.N.I. sensor is a disposable, pre-gelled 3-lead array. Place the large and the small electrodes on the patient in order to detect a cardiac signal of good amplitude. For example, the large electrode can be put on the right shoulder and the small electrode at the apex of the heart as shown in the diagram on the pouch label.

The sensor is self-adhesive and contains bio-compatible materials where the array comes into contact with the skin.

The sensor is designed for use in the Adult population and for the population above the age of 2 years although supervision is required.

The A.N.I. sensor is to be used with the A.N.I. Pain Monitoring System and may be worn for a continuous period of 24 hours in a hospital environment.

Cautions:

Rx only – US Federal law restricts this device to sale by or on the order of a physician.

A.N.I. sensor may damage the skin if removed carelessly.

The A.N.I sensor should be applied only to intact, clean skin (e.g. not over open wounds, lesions, infected, or inflamed areas).

The A.N.I sensor should only be used by or in consultation with a healthcare provider familiar with their proper placement and use.

Ensure the skin area is clean, dry and hair free before sensor application. Use alcohol wipes for cleaning the application area.

Apply the A.N.I. sensor within 1 hr of removal from the inner packaging.

Replace the A.N.I. sensor if it no longer sticks firmly to the skin.

Reddening of the skin and slight irritation may occur. This is normal.

Dispose of the A.N.I sensor properly after use under hospital guidelines.

The sensor is not waterproof.

Store in a cool dry place.

The A.N.I sensor is not made with natural rubber latex.

Precautions:

During surgical procedures, the A.N.I. sensor should be placed as far as possible from any electro-surgical area to minimise unwanted RF current flow.

Do not wear the A.N.I sensor during an MRI or X-Ray

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